



1064nm/ 532nm Q-switched Nd: YAG Laser

Genesis-Q utilizes two wavelengths 1064nm/ 532nm

- and deliver laser power in a very short period of time
- as triple pulsed & minimize treatment downtime.
 Auto Calibration & Self Restoration.
- Energy is delivered efficiently to maximize effect while minimizing damage to surrounding tissues.
- Strong Laser toning function by high energy & efficient delivery system.
 - Dual Laser 1064nm/ 532nm fractional function for
- various operation.
- Application for an effective power supply & cooling system.
 Design for a stable resonator to minimize an external environment effect.

Applications

Melasma
Epidermal Nevi
Tattoo Removal
Active Acne
Non-Ablative resurfacing
Pigminted & Vascular
lesons and many more..


Options

For user convenience, the system recognizes the hand piece among three different types, One Zoom Hand-piece & two Fractional Hand-piece, and its spot size automatically.

Two different types of Hand-piece tips are also provided for more comfortable operation.



- Deeper Penetration depth
 Multi pulse laser beam penetrates
 to dermis
- Patients need less time to remove tattoos, Pigmented lesions and vascular lesions
- Painless treatment
 Small peak of multi pulse results in less
 pain for patients.
- No need of anesthesia







1064nm Fractional Hand-piece

9x9mm2 / 9x9 dots

532nm Fractional Hand-piece















System Specification

General

Dimension (w x D x H) 420(W) x 805(D) x 1620(H) mm 90 KG Weight







Nd: YAG Laser

Genesis-Q

